



ECOCARE SERVICES LLP
My SingPost Box 880171
Singapore 919191
Tel/Fax: 65 63103163
Mobile: 65 96857026
Email:ecocare_services@yahoo.com.sg
UEN No:T07LL0059J
GST Reg No:M90357607P

A-318 Disinfectant

DESCRIPTIONS: **A-318 Disinfectant** is a **Hospital-grade** disinfectant, deodorant, fungicide, and virucide that reduces and controls airborne bacteria and virus particles. It is effective against a broad spectrum of gram negative and gram positive organisms that helps to prevent cross-contamination of germs and bacteria

APPLICATIONS: **A-318 Disinfectant** is widely used in hospitals, clinics, schools, institutions, hotels, nursing homes, restaurant, offices, commercial areas etc. It can also be used to disinfect and deodorize all hard surfaces including table-tops, chairs, toys, toilets, floors, walls, air-conditioning etc.

ADVANTAGES:

- **CONTROLS SPREAD OF VIRUSES**
A-318 Disinfectant not only aids in controlling and reducing spread of virus infections but also reduces the hazard of cross-contaminations.
- **HIGHLY EFFECTIVE DISINFECTANT AND SANITIZER**
It kills germs on surfaces with its powerful actions that work on frequently handled objects which may spread disease-causing particles.
- **ELIMINATE ODOUR**
Beside sanitizing and disinfecting, **A-318 Disinfectant** also eliminates odour, leaving areas clean and gives a lingering scent after application.
- **MILDEW CONTROL**
A-318 Disinfectant helps to prevent growth of mold and mildew on surfaces.
- **ECONOMICAL**
It can be diluted with water for cost saving.

EFFECTIVE:

Its powerful disinfecting action works on frequently handled objects which may spread disease-causing organism. If used as directed, it also disinfects pre-cleaned, hard non-porous surfaces against Staphylococcus Aureus, Salmonella Choleraesuis and Pseudomonas Aeruginosa, thus meeting requirements for hospital use. **A-318 Disinfectant** is also effective against other common viruses such as: Influenza, Severe Acute Respiratory Syndrome (SARS), MERS, Coronavirus, HFMD and surrogates, etc..

DIRECTIONS:

- **A-318** may be applied by mopping, brushing, soaking, circulating and mist spraying. For disinfecting purpose, it is recommended to dilute 1 part of **A-318** to 25 parts of water.
- For sanitizing purpose, dilution rate is 1:40 - 80 parts of water. Brush the entire surface to remove soil and films, and rinse away with water.
- For Hard Surfaces e.g. table tops, chairs sinks, cabinets etc: Recommend weekly treatment for sanitizing and disinfecting. Spray over surfaces, let contact stay for 3-5 minutes before wiping off.

PROPERTIES:

Appearance: Purple red liquid.
Odour: Floral scent.
Active Ingredient: Quaternary ammonium compound with biodegradable nonionic surfactant